

## Portable Floor Model Surgery Table

The Shank's® Portable Floor Model Surgery Table was introduced in 1978. This table was the first mobile large animal surgery table that incorporated a scissor-type lift design. The Table can be placed in a "stationary" position by retracting the casters, and is also capable of 2-way tilt functions.

### Technical Specs

- 36" x 72" top (91cm x 183cm)
- 6" swivel casters (with swivel lock capability)
- Self-contained power unit (available 110V/60 cycle or other power configurations for international use)
- Hand control hydraulic valve (3 functions)
- 2-way hydraulic tilt capability
- Versatile removable supports for both dorsal and lateral positioning of animal
- Lowered height 14" (36cm)
- Raised height 48" (122cm)
- Hot dip galvanized coating/enamel topcoat
- Stainless steel fasteners & motor shield
- Rated for 2,500 lbs. (1,134 kg)

### Standard Accessories

- 3 flat extensions (22" x 26"/56cm x 66cm)
- Offset bracket for offset placement of flat extension for head

### Optional Accessories/Features

- 36" x 84" top (91cm x 213cm)
- Removable cutouts (on table top)
- 1-1/2" (3.8cm) thick removable Ensolite® pad
- V-head extension
- Small shoulder extension
- Leg poles, dorsal padded props & tie loops
- SS intestinal tray
- SS intestinal tray with free standing base
- Lateral arthroscopic leg bracket
- Simple lateral bracket
- Dorsal arthroscopic leg bracket (pair)

- SS IV Pole (with or without breathing circuit)
- Round tube cross pipe
- T-handle
- P/S storage cart
- Foam positioning pads
- Lifting hobbles

### Hydraulic System and Operation

The Table is user friendly and merely requires access to power for operation. Once the table is connected to power, an on/off switch turns on the hydraulic system. The Table is raised and lowered by manipulation of the hand control lever mounted on the table. The table can be stopped at any position during raising or lowering. The remaining two handles allow for independent retraction of the casters on the head and tail end of the table.

The table is completely mobile within a clinic setting via 6" (15cm) casters.

All optional accessories are easily connected to the table by inserting them into any of the table ports located around the perimeter of the table top. The optional accessories allow for both dorsal and lateral support. Accessories are available to be used in any and all surgery procedures commonly performed.

