

CT Platform Table

The Shank's® CT Platform Table ("Table") was introduced in 2019 and is designed for support of standing large animals in conjunction with CT scanning. The Table is designed to be placed in a recessed pit adjacent to the CT machine.

Technical Specifications

- 60-1/4" x 96-7/8" Table top with non-slip rubber surfacing
- Table includes an affixed Stock Frame
- Table can raise 26-1/2" above surrounding floor
- Table can lower 31" into pit
- Self-contained power unit (available 110V/60 cycle or other power configurations for international use)
- Hot dip galvanize coating/enamel topcoat
- Rated for 2,500 lbs. (1,134 kg)

Stock Frame/Structure

- 2" pipe frame
- Rear locking gate entry
- Swinging side rails
- Front slide gate
- Adjustable breast bar
- Adjustable back up bar with pad

Hydraulic System and Operation

This Table requires on site hydraulic system installation and anchoring inside the pit to complete the installation. Once installed, the Table is raised and lowered by a hand held push button pendant with UP and DOWN functions. The table can be stopped at any position during raising or lowering.

